




Skin and Wound Product Information Sheet

Alldress	
Classification	Cover Dressing: Composite
Key Points	<ul style="list-style-type: none"> • A multilayer, absorbent, vapour-permeable dressing with border; border is coated with a water-based solvent-free polyacrylate adhesive • Laminated covering of vapour permeable polyurethane film that acts as a barrier to fluid strike-through, bacteria and viruses • Moisture retentive • Showerproof • Smooth surface to reduce friction
Indications	<ul style="list-style-type: none"> • For wounds with small to moderate amounts of exudate
Precautions	<ul style="list-style-type: none"> • Use with caution on fragile skin
Contraindications	<ul style="list-style-type: none"> • Do not use on untreated clinically-infected wounds
Formats & Sizes	<ul style="list-style-type: none"> • Adhesive border dressing <ul style="list-style-type: none"> ▪ 10 x 10 cm ▪ 15 x 15 cm ▪ 15 x 20 cm <div style="text-align: right; margin-top: 10px;">  </div>
Application Directions	
Cleanse the wound and the peri-wound skin with sterile normal saline or agency approved wound cleanser; dry peri-wound skin.	Reduces wound debris and allows for adhesion of dressing or tape.
If required, apply skin barrier to peri-wound skin.	To protect the peri-wound skin from maceration and to improve the adhesion of the dressing or tape.
To Apply	
Apply directly on the wound as a primary dressing or use in combination with another wound product as a secondary cover dressing.	
Corners may be slit to accommodate body parts eg elbows, heels	
Smooth adhesive border to the skin.	To ensure good adhesion to the skin.
Do not stretch the dressing.	To prevent client discomfort with movement of the affected body part.
To Remove	
Gently lift the border of the dressing.	To minimize trauma to the peri-wound skin.
Frequency of Dressing Change	
Will depend upon the amount of exudate. May be left in place 1 to 5 days.	
Expected Outcome	
Exudate is managed with no peri-wound skin maceration.	
For further information, please contact your Wound Clinician.	