




Skin and Wound Product Information Sheet

Exu-dry	
Classification	Contact Layer: Anti-shear
Key Points	<ul style="list-style-type: none"> • Sterile absorbent anti-shear dressing
Indications	<ul style="list-style-type: none"> • For wounds with moderate to large amounts of exudates • For painful wounds, graft and donor sites
Precautions	<ul style="list-style-type: none"> • N/A
Contraindications	<ul style="list-style-type: none"> • N/A
Formats & Sizes	<ul style="list-style-type: none"> • Sheet <ul style="list-style-type: none"> ○ 10 x 10cm ○ 15 x 23cm ○ 23 x 38cm ○ 38 x 46cm ○ 38 x 61cm
	
Application Directions	Rationale
Cleanse/irrigate wound with sterile normal saline or agency approved wound cleanser; dry peri-wound skin.	Reduces wound debris and allows for adhesion of dressing or tape.
If required, apply skin barrier to peri-wound skin.	To protect the peri-wound skin from maceration and to improve the adhesion of the dressing or tape.
To Apply	
Choose appropriate size of dressing. Do not cut dressing.	
Apply directly to wound; secure in place with tape, kling or mesh stocking.	
To Remove	
Gently remove the dressing from the wound bed.	To avoid trauma to the wound bed.
If adhered, soak with sterile normal saline or remove in shower /tub. The dressing will soften and become heavy as it absorbs moisture; remove gently.	
Frequency of Dressing Change	
Will depend upon the amount of exudates. Dressing may be left on wound for up to 7 days.	
Expected Outcome	
Exudate is managed with no peri-wound skin maceration.	
For further information, please contact your Wound Clinician.	