









Skin and Wound Product Information Sheet

| Mepilex Border for Pressure Injury Prevention | |
|--|---|
| Classification | Prophylactic Dressing |
| Key Points | <ul style="list-style-type: none"> Used to enhance, but not replace, routine pressure injury prevention strategies for the prevention of sacral/coccyx and heel pressure injuries. Highly conformable 5-layer foam dressing with Deep Defence™ technology that redistributes pressure/shear forces, reduces friction, and balances microclimate. Soft and conformable waterproof/showerproof foam dressing with silicone adhesive layer for atraumatic dressing removal. |
| Indications | <ul style="list-style-type: none"> For those patients, clients or residents at risk for developing a pressure injury related to pressure, friction/shear on the sacral-coccyx or heel area. |
| Precautions | <ul style="list-style-type: none"> Does not replace the use of other pressure injury prevention strategies (i.e. pressure risk assessment, regular positioning, appropriate pressure redistribution and support strategies). Consult with Physician/NP/Wound Clinician prior to using foam dressings (of any kind) on ischemic lower legs/feet Dressing should not be cut |
| Contraindications | <ul style="list-style-type: none"> Pre-existing pressure injury on the sacral-coccyx or heel, including Pressure Injury Stage 1 Trauma or burn to sacrum or coccyx. Do not use with skin barriers /skin sealants or cleansing wipes containing dimethicone/silicone, emollients etc. as these reduce the effectiveness of the adhesive properties of the dressing |
| Formats & Sizes | <ul style="list-style-type: none"> Sacrum <ul style="list-style-type: none"> ○ 16 x 10cm ○ 22 x 25cm Heel <ul style="list-style-type: none"> ○ 18 x 24cm <div style="text-align: center;">  </div> |
| Application Directions | |
| To Apply to the Sacral-coccyx Area | |
| <p>Cleanse the sacral area with pH-balanced skin cleanser or warm water. Gently pat the skin dry. Do not use emollient or dimethicone/silicone ointments/barrier wipes or skin sealants in area where dressing will be applied.</p> <p>Have a colleague hold the buttocks apart. Remove the dressing's center release film and apply dressing into the gluteal cleft and then sacral area. Remove the right-side release film and gently smooth this side of the dressing into place. Repeat with the left side. Run the side of a hand along the gluteal cleft to secure placement.</p> | <p>Emollients, dimethicone/silicone (and other skin preparations) can reduce the adhesive properties of the silicone dressing.</p> <p>It is important that the dressing 'fits' into the upper aspect of the gluteal cleft to ensure that the dressing is properly secured against incontinence episodes.</p> |
| <div style="display: flex; justify-content: space-around;">   </div> | |
| <p>On the dressing, print "P" for Preventative dressing and add the date that the dressing was applied.</p> | <p>To communicate with other staff the purpose of the dressing and when it needs to be changed</p> |
| To Apply to the Heel Area | |
| <p>Apply the adherent part of the dressing marked 'A' to the posterior heel/Achilles tendon areas, positioning the narrowest part of the dressing at the base of the heel. Do not stretch.</p> <p>Remove the backing from one of the ankle flaps. Apply and smooth. Repeat with the other side.</p> | <div style="display: flex; justify-content: space-around;">   </div> |



Skin and Wound Product Information Sheet

| Application Directions | Rationale | |
|--|--|--|
| <p>Gently apply the adherent part of the dressing marked 'B' under the plantar surface of the foot. Do not stretch.</p> <p>Remove the backing from one of the flaps with tabs. Apply and smooth border. Repeat with the other side.</p> <p>Press and smooth the dressing to ensure the entire dressing is in contact with the skin.</p> |  | |
| <p>On the dressing, print "P" for Preventative dressing and add the date that the dressing was applied.</p> | <p>To communicate with other staff the purpose of the dressing and when it needs to be changed.</p> | |
| Daily Care | | |
| <p>As part of the evaluation of client's specific pressure injury prevention strategies, once a day, peel dressing back and assess the skin. Reapply existing dressing ensuring the border of the dressing is smooth with no wrinkles. Document assessment findings.</p> <p>If patient is incontinent and top of dressing is soiled, gently wipe off.</p> | <p>Wrinkles in the dressing are due to shear forces being applied to the dressing; if possible, remove these concerns e.g. lower the head of the bed.</p> <p>Dressing is waterproof and will not allow urine or feces to soak into the dressing.</p> | |
| <p>If dressing does not stay intact for longer than 24 hours due to incontinence, discontinue the dressing and use barrier cream or alternative skin management.</p> <p>If a pressure ulcer develops within the area of the dressing, discontinue the prevention/protection dressing and initiate appropriate wound management. Inform OT, PT and/or Wound Clinician of the pressure ulcer occurrence.</p> | <p>If dressing does not stay in place it is not a cost-effective prevention strategy.</p> <p>Other pressure ulcer interventions will need to be considered.</p> | |
| To Remove | | |
| <p>Gently lift the border and use one hand to stabilize the skin.</p> | <p>To minimize trauma to skin.</p> | |
| Frequency of Dressing Change | | |
| <p>May be left in place for up to 7 days. Change the dressing should it lose its adherence e.g. edges roll, border does not 'stick', or become soiled.</p> <p>Replace the prevention dressing as long as patient meets selection criteria above.</p> | <p>As the patient's level of pressure ulcer risk improves, the dressing may no longer be required.</p> | |
| Expected Outcome | | |
| <p>Pressure injury does not develop.</p> | | |
| <p>For further information, please contact your Wound Clinician.</p> | | |