




## Skin and Wound Product Information Sheet

<b>Mepitel One</b>	
<b>Classification</b>	<b>Layer: Silicone</b>
<b>Key Points</b>	<ul style="list-style-type: none"> <li>One sided wound contact layer which conforms to body contours and does not leave a residue or dries out</li> <li>Allows for cover dressing to be changed as needed without disruption of the wound bed</li> <li>Does not stick to gloves or scissors</li> <li>Low potential for skin irritation/allergy</li> <li>Can be left in place for up to 14 days if there are no signs of infection</li> </ul>
<b>Indications</b>	<ul style="list-style-type: none"> <li>For protection of a healing wound bed from secondary dressing</li> <li>For painful wounds</li> <li>May used in conjunction with Negative Pressure Wound Therapy (NPWT)</li> </ul>
<b>Precautions</b>	<ul style="list-style-type: none"> <li>If used for Epidermolysis Bullosa, remove carefully as adhesion level is slightly higher than regular Mepitel</li> <li>If used for burns treated with meshed grafts or after facial resurfacing, imprinting of contact layer can occur; avoid unnecessary pressure upon the dressing and reposition the contact layer at least every second day</li> <li>If used for the fixation of skin grafts, the dressing should not be changed before the fifth day post application or as per surgeon's order</li> </ul>
<b>Contraindications</b>	<ul style="list-style-type: none"> <li>Do not use with skin sealant/barrier</li> </ul>
<b>Formats &amp; Sizes</b>	<ul style="list-style-type: none"> <li>Sheet                             <ul style="list-style-type: none"> <li>5 x 7.5 cm</li> <li>7.5 x 10 cm</li> <li>10 x 18 cm</li> <li>17 x 25 cm</li> </ul> </li> </ul>
	
Application Directions	Rationale
Cleanse/irrigate wound with sterile normal saline or agency approved wound cleanser; dry peri-wound skin.	Reduces wound debris and allows for adhesion of dressing or tape.
<b>To Apply</b>	
Remove the release film.	To prepare dressing for application.
Apply Mepitel One, tacky side to the wound.	To ensure the correct side is to the wound.
Apply appropriate cover dressing to maintain a moisture-balanced wound environment.	The choice of cover dressing is depended upon the amount of exudate expected.
<b>To Remove</b>	
Gently lift layer off of the wound bed	To avoid trauma the wound bed.
<b>Frequency of Dressing Change</b>	
Will depend upon the amount of exudate. Mepitel One may be left on the wound for up to 14 days.  Reposition the contact layer at least every second day when used for burns treated with meshed grafts or after facial resurfacing.  If using as contact layer under NWPT, then apply new contact layer with each dressing change.	Dressing is transparent enough to allow for assessment of wound.
<b>Expected Outcome</b>	
Contact layer does not adhere to wound bed.	
<b>For further information, please contact your Wound Clinician.</b>	