


## Skin and Wound Product Information Sheet

<b>TIELLE NON-ADHESIVE</b>	
<b>Classification</b>	<b>Cover Dressing: Hydropolymer Foam</b>
<b>Key Points</b>	<ul style="list-style-type: none"> <li>• A non-border hydropolymer dressing composed of three layers: a wound contact layer, a superabsorbent foam layer and a polyurethane film backing</li> <li>• Flexible</li> <li>• Can be cut to shape</li> <li>• Needs to be secured</li> <li>• May be used as a primary or secondary dressing</li> </ul>
<b>Indications</b>	<ul style="list-style-type: none"> <li>• For wounds which have small to large/copious exudate</li> <li>• May be used under compression wraps</li> </ul>
<b>Precautions</b>	<ul style="list-style-type: none"> <li>• May be used on wounds with S&amp;S of infection when the infection is being treated</li> </ul>
<b>Contraindications</b>	<ul style="list-style-type: none"> <li>• Do not use on lesions with active vasculitis</li> <li>• Do not use on full thickness burns</li> </ul>
<b>Formats &amp; Sizes</b>	<ul style="list-style-type: none"> <li>• Dressing                             <ul style="list-style-type: none"> <li>▪ 5 x 5cm</li> <li>▪ 10 x 10cm</li> <li>▪ 15 x 15cm</li> <li>▪ 17 x 17cm</li> <li>▪ 21 x 22cm</li> </ul> </li> </ul> <div style="text-align: right; margin-top: 10px;">  </div>
<b>Application Directions</b>	
Choose a dressing which will ensure at least a 1cm overlap onto the peri-wound skin.	<b>Rationale</b>
Cleanse the wound/peri-wound skin with at least 100ml of sterile normal saline or agency approved wound cleanser; dry peri-wound skin.	Reduces wound debris and allows for adhesion of dressing or tape.
If required, apply skin barrier to peri-wound skin.	To protect the peri-wound skin from maceration and to improve the adhesion of the dressing or tape.
<b>To Apply</b>	
If needed, cut dressing to fit wound area.  Partially peel off the backing paper. Position the dressing on the wound and ensure that there is 1cm overlap on the peri-wound skin.  Remove remaining backing paper, smooth into place and secure with tape or Kling/tape.	To hold dressing in place.
<b>To Remove</b>	
Gently remove tape or Kling/tape and lift one corner of the dressing and remove dressing.	To minimize trauma to the peri-wound skin.
<b>Frequency of Dressing Change</b>	
Change the dressing when the exudate is within 1cm of the edge of the foam. Do not allow exudate to accumulate under the backing.  Depending upon the amount of exudate, the dressing may be left in place for up to 7 days.	
<b>Expected Outcome</b>	
Exudate is managed with no peri-wound skin maceration.	
<b>For further information, please contact your Wound Clinician.</b>	