




Skin and Wound Product Information Sheet

UrgoTul	
Classification	Contact Layer
Key Points	<ul style="list-style-type: none"> A conformable, flexible, non-occlusive, and non-adherent contact layer comprised of a lipido-colloid matrix impregnated onto a fine polyester mesh. The lipido-colloid matrix forms a gel in the presence of wound exudate to provide a moist interface between the wound and the dressing. Helps reduce or eliminate pain at dressing changes and minimizes damage to newly formed tissue. Has a gentle tack to help hold dressing in place. Requires a cover dressing. May be left in place for up to 7 days.
Indications	<ul style="list-style-type: none"> Partial and full thickness wounds with low to moderate amounts of exudate. May be used as an interface under negative pressure wound therapy dressing. May be used on mobile joints.
Precautions	<ul style="list-style-type: none"> N/A
Contraindications	<ul style="list-style-type: none"> Those with sensitivity or allergy to any dressing component.
Formats & Sizes	<ul style="list-style-type: none"> Sheet <ul style="list-style-type: none"> 5 x 5 cm 10 x 12.5 cm 15 x 20 cm
	
Application Directions	Rationale
Cleanse/irrigate wound with sterile normal saline or agency approved wound cleanser.	Reduces wound debris and allows for a adhesion of dressing or tape.
Clean and dry periwound skin.	
To Apply	
Choose a size that covers the entire wound. The dressing may overlap onto periwound skin or it may be cut to size. When cutting the dressing, ensure that the protective backing in place.	Overlap onto skin assists with adhesion and wear time.
Use latex free gloves or moisten latex gloves with sterile saline prior to handling dressing.	Dressing tends to stick to latex gloves.
Remove the protective clear film from both sides of the dressing.	
Apply the dressing to cover the entire wound bed. The dressing may overlap onto healthy skin.	
Apply appropriate cover dressing to maintain a moisture-balanced wound environment.	The choice of cover dressing is dependent upon the amount of exudate expected.
To Remove	
Using forceps, gently remove dressing from the wound.	To avoid trauma to the wound.
Frequency of Dressing Change	
Change the contact layer at least every 7 days. Change cover dressing as needed.	
Expected Outcome	
A moist wound environment is maintained.	
No pain present upon removal of dressing.	
For further information, please contact your Wound Clinician.	