





Product Information Sheet

Biatain Alginate Backorder Substitution for Kaltostat	
Classification	Wound Fillers: Calcium Alginates
Key Points	<ul style="list-style-type: none"> Biatain Alginate is being used as a PHSA Supply Chain approved substitution for Kaltostat sizes 5x5 cm and 7.5x12 cm. An absorbent, non-woven pad comprised of 85% calcium alginate and 15% carboxymethylcellulose that forms a gel as wound exudate is absorbed. Requires a secondary dressing.
Indications Precautions Contraindications	<ul style="list-style-type: none"> Refer to Kaltostat Product Information Sheet. Do not use on those with a known sensitivity or allergy to calcium alginate or carboxymethylcellulose.
Formats & Sizes	<ul style="list-style-type: none"> Sheet: <ul style="list-style-type: none"> 5x5 cm 10x10 cm 15x15 cm Rope: <ul style="list-style-type: none"> 3x44cm <div style="display: flex; justify-content: center; align-items: center; gap: 20px;">   </div>
Application Directions	<ul style="list-style-type: none"> Refer to Kaltostat Product Information Sheet.
Dressing Change Frequency	<ul style="list-style-type: none"> Refer to Kaltostat Product Information Sheet.
For further information please contact NSWOC/Wound Clinician	