Created by the British Columbia Provincial Professional Practice Stream – Wound Ostomy Continence







•









Providence Health Care

Medigrip LF Lower Limb Application for Compression			
Limb Measurement	Compression**/Size of Bandage		
Calf Circumference* Minimum - Maximum Centimeters	Low 5-10mmHg	Medium 10-20mmHg	High 20-30mmHg
13.0 - 14.0cm	А	-	-
14.0 - 15.0cm	В	А	-
15.0 - 18.0cm	С	В	А
18.0 - 20.0cm	D	В	А
20.0 - 23.0cm	E	С	В
23.0 - 25.0cm	F	D	В
25.0 - 30.0cm	F	E	С
30.0 - 38.0cm	G	F	D
38.0 - 46.0cm	G	F	E
46.0 - 58.0cm	J	G	F
58.0 - 71.0cm	J	G	-
71.0 - 97.0cm	К	J	-
*Measure widest part of the calf. ** Compression mmHg achieved with double layer of stocking.			