




## Skin and Wound Product Information Sheet

<b>TIELLE Plus</b>	
<b>Classification</b>	<b>Cover Dressing: Hydropolymer Foam</b>
<b>Key Points</b>	<ul style="list-style-type: none"> <li>• An adhesive foam dressing composed of three layers: a hydropolymer foam, a superabsorbent wicking layer and polyurethane backing.</li> <li>• Has an adhesive border</li> <li>• Flexible</li> <li>• May be used as either a primary or secondary dressing.</li> </ul>
<b>Indications</b>	<ul style="list-style-type: none"> <li>• For wounds which have moderate to large/copious amount of exudate</li> <li>• May be used under compression wraps</li> </ul>
<b>Precautions</b>	<ul style="list-style-type: none"> <li>• May be used on wounds with S&amp;S of infection when the infection is being treated.</li> </ul>
<b>Contraindications</b>	<ul style="list-style-type: none"> <li>• Do not use on full thickness burns</li> <li>• Do not use on lesions with active vasculitis</li> </ul>
<b>Formats &amp; Sizes</b>	<ul style="list-style-type: none"> <li>• Dressing size/Pad Size                             <ul style="list-style-type: none"> <li>▪ 11 x 11cm/7.2 x 7.2cm</li> <li>▪ 15 x 15cm/11 x 11cm</li> <li>▪ 15 x 20cm/11 x 16cm</li> <li>▪ 15 x 15cm (sacral)/10 x 7.2cm</li> <li>▪ 20 x 26.5cm (heel)/14.5 x14.5cm</li> </ul> </li> </ul> 
<b>Application Directions</b>	
Choose a dressing which has a pad size that will overlap onto the peri-wound skin by at least 1cm.	
Cleanse the wound/peri-wound skin at least 100ml of sterile normal saline or agency approved wound cleanser; dry peri-wound skin.	Reduces wound debris and allows for adhesion of dressing or tape.
If required, apply skin barrier to peri-wound skin.	To protect the peri-wound skin from maceration and to improve the adhesion of the dressing or tape.
<b>To Apply</b>	
Partially peel back the side backing papers.	Prepares dressing for application.
Position the pad over the wound and ensure that there is a 1cm overlap on the peri-wound skin. Smooth in place.	Ensures dressing is aligned with wound so that exudate is effectively absorbed.
Peel off the backing papers and smooth adhesive border onto peri-wound skin.	Provides adequate seal.
Adding petroleum jelly over the dressing borders will add additional water repellency if required.	
<b>To Remove</b>	
Gently lift one corner of the dressing and gently remove dressing.	To minimize trauma to the peri-wound skin.
For fragile peri-wound skin, lift one corner of the dressing and apply a normal saline soaked Q-tip to the underside of the dressing as it is peeled back.	To break the adhesive seal.
<b>Frequency of Dressing Change</b>	
Change the dressing when the exudate is present at the edge of the foam pad. Do not allow exudate to accumulate under the backing.	
Depending upon the amount of exudate, the dressing may be left in place for up to 7 days.	
<b>Expected Outcome</b>	
Exudate is managed with no peri-wound skin maceration.	
<b>For further information, please contact your Wound Clinician.</b>	