




Product Information Sheet

Ultrisorbs Under Pad	
Classification	Continence Containment: Linen Protector
British Columbia Practice	<ul style="list-style-type: none"> Do not use under pad as a patient positioning or lifting device as per Health Authority's Occupational Health & Safety guidelines.
Key Points	<ul style="list-style-type: none"> Disposable under pad designed with a: <ul style="list-style-type: none"> White breathable top sheet. Middle layer of super absorbent polymer (SAP) material which can manage light to heavy amounts of fluid with moisture being absorbed up to 5cms from the edges of the pad. Also provides odour management. Back sheet that is air-permeable and leak-proof until the pad is saturated. The air-permeable feature allows the under pad to be used on a support surface with microclimate features (e.g., low air-loss mattress or microclimate mattress cover, i.e., a Skin IQ). Do not layer under pads. Do not use as a procedure pad.
Indications	<ul style="list-style-type: none"> For the client who is not using a continence containment product, (e.g., a pull-up brief) and <ul style="list-style-type: none"> Is undergoing toilet retraining and is at risk for a urinary and/or fecal incontinent episode. Has stress or urge incontinence and is at risk for small urinary incontinent episode. Do not use for the client who is using a or a continence device, (e.g., a urinary catheter) to minimize the number of layers beneath the client.
Precautions	<ul style="list-style-type: none"> N/A.
Contraindications	<ul style="list-style-type: none"> Sensitivity or allergy to any of the components of the pad.
Formats & Sizes	<ul style="list-style-type: none"> Pad Size: <ul style="list-style-type: none"> 57.5 x 90 cm 75.0 x 90 cm <div style="text-align: right;">  </div>

Directions	Rationale / Key Points
Selection	
Select size of pad needed for bed or chair.	
Preparation	
Ensure the bed or chair's bottom sheet is not wrinkled.	
Application	
Lay the pad down with the white side up. Smooth out to remove any wrinkles. Use only one pad. Do not fold or tuck pad into the perineal area.	The more layers of linen (e.g., bottom sheet, linen protector, positioning slings/devices, etc.), between the client and the support surface, the higher the risk of the client developing a pressure injury. Under pad must lie flat to collect and 'gel' the fluid.
Removal	
Do not pull the pad out from under the client.	Pulling the pad may cause the pad to tear and/or a Medical Device Related Skin Damage injury.
Frequency of Dressing Change	
Change when soiled. Product is disposable; do not launder.	Should the pad reach its saturation point, the pad no longer pulls fluid into its core. The client's skin is now in



Product Information Sheet

Directions	Rationale / Key Points
	contact with the fluid, (e.g., urine) leading to the risk for Incontinence Associated Dermatitis (IAD).
Expected Outcomes	
<p>The under pad will protect the linen from incontinence episodes.</p> <p>The client will not develop Incontinence Associated Dermatitis (IAD).</p> <p>Product performs as expected.</p>	<p>If product does not perform as expected, notify NSWOC/Wound Clinician and then consider submitting a Supply Chain Product Concern Form.</p>
For further information please contact NSWOC/Wound Clinician	