




Product Information Sheet

CURAD Non-Adherent Oil Emulsion Dressing Backorder Substitution for Adaptic	
Classification	Category: Subcategory
Key Points	<ul style="list-style-type: none"> • CURAD Non-Adherent Oil Emulsion Dressings are being used as a PHSA Supply Chain approved substitution for Adaptic sizes 12.7 x 22.9 cm and 7.6 x 152.4 cm • A knitted, high porosity cellulose acetate dressing impregnated with white petrolatum. • Primary dressing requiring a secondary dressing.
Indications Precautions Contraindications	<ul style="list-style-type: none"> • Refer to Adaptic Product Information Sheet
Formats & Sizes	<ul style="list-style-type: none"> • Sheet: <ul style="list-style-type: none"> ○ 7.6 x 7.6 cm ○ 7.6 x 20.3 cm (1 per pkg) ○ 7.6 x 20.3 cm (3 per pkg) ○ 7.6 x 40.6 cm 
Application Directions	<ul style="list-style-type: none"> • Refer to Adaptic Product Information Sheet
Dressing Change Frequency	<ul style="list-style-type: none"> • Daily or as instructed by Physician, NP or NSWOC/Wound Clinician
For further information please contact NSWOC/Wound Clinician	